

Gdańsk, 23. December 2020

To whom it may concern,

On behalf of the GeneMe, taking into account the latest information on the mutation of the S-gene and its influence on the results of RT-PCR assays, we hereby present our research and development report, the objective of which was to determine if the recently observed mutations in the SARS-CoV-2 affect the loss of specificity of the FRANKD RT-LAMP test. This research and development report summarizes our findings regarding the influence of S gene mutations in SARS-CoV-2 on FRANKD primers hybridization. For this purpose, bioinformatic analysis of the mutated sequences of the coronavirus S gene was performed to assess if the point mutation or deletions lay in the hybridization region of FRANKD RT-LAMP primers.

RT-LAMP technologies are less prone to lose specificity and sensitivity in comparison to RT-PCR, because of (i) amplification principle and reaction efficacy (constant amplification, no doubling DNA in each cycle), (ii) more than two primers in the reaction (more targets), (iii) RT-LAMP polymerase features (namely strand display activity of the polymerase).


I can definitively state that FRANKD RT-LAMP assay's ability to detect SARS-CoV-2 remains at the highest level regardless of this new mutation.


GeneMe is working with the head of RT-LAMP in the UK Dr Veronica Fowler and Dr Lesley Scott, the Head of the National Laboratory in South Africa, to ensure a direct line into the latest research from these two countries.

Yours faithfully,

Kasjan Szemiako, Research and Development Director, GeneMe





| | |
|--|--|
|  <p>G E N E M E</p> | <p>Research & Development</p> |
| <p>RESEARCH & DEVELOPMENT REPORT</p> | |
| <p>NAME AND SURNAME</p> | <p>Head of RD Department at GeneMe, Marta Skwarecka, PhD</p> |
| <p>ANALYZED IVD PRODUCT</p> | <p>FRANKD</p> |
| <p>DATE OF REPORT PREPARATION</p> | <p>22.12.2020</p> |
| <p>TITLE OF THE REPORT</p> | <p>The analysis of mutations of S gene in SARS-CoV-2 and their influence on FRANKD primers hybridization</p> |
| <p>1. Methods</p> <p>The analysis aimed to show if the recently observed mutations in the SARS-CoV-2 virus genome [1-5] affect the FRANKD test's specificity and assess the potential reduction of efficient hybridization of the 6 primers used in FRANKD.</p> <p>Bioinformatic analysis of the mutated sequences of the coronavirus S gene was performed, with particular emphasis on nucleotide primers sequences.</p> <p>The literature review indicated mutations within the S gene: deletion 69-70, deletion 144, N501Y, N439K, D614G, A570D, P681H, T716I, S982A, D1118H (date of accession: 23.12.2020, GeneBank NCBI).</p> <p>The mutations were analyzed if they lay in the region of FRANKD primers hybridization.</p> | |
| <p>2. Results</p> <p>Below is presented the nucleotide sequence of the wild-type SARS-CoV-2 S gene with the point mutations marked in red (date of analysis 22 December 2020 – current stage of described mutations in the scientific literature – NCBI GeneBank database). NCBI GeneBank database collects data from all published research papers and GISAID EpiCov database and therefore ensures consistency.</p> | |

ATGTTTGTTCCTTGTTCCTTATTGCCACTAGTCTCTAGTCAGTGTGTTAATCTTACAACCAGAACTCAATT
ACCCCTGCATACACTAATTCTTTCACACGTGGTGTTCCTTACCCTGACAAAAGTTTTCAGATCCTCAGTTT
TACATTCAACTCAGGACTTGTTCCTTACCTTTCTTTTCCAATGTTACTTGGTTCATGCTATA**CATGTC**TC
TGGGACCAATGGTACTAAGAGGTTTGATAACCCTGTCTACCATTTAATGATGGTGTTCATTTTGCTTCCA
CTGAGAAGTCTAACATAATAAGAGGCTGGATTTTGGTACTACTTTAGATTTCGAAGACCCAGTCCCTACTT
ATTGTTAATAACGCTACTAATGTTGTTATTAAGTCTGTGAATTTCAATTTTGAATGATCCATTTTGGG
TGTT**TAT**TACCACAAAAACAACAAAGTTGGATGGAAAGTGAGTTCAGAGTTTATTCTAGTGCGAATAAT
TGCACCTTTGAATATGTCTCTCAGCCTTTTCTTATGGACCTTGAAGGAAAAACAGGGTAATTTCAAAAAATCT
TAGGGAATTTGTGTTAAGAATATTGATGGTTATTTTAAAAATATATTCTAAGCACACGCCTATTAATTTAG
TGCGTGATCTCCCTCAGGGTTTTTTCGGCTTTAGAACCATTGGTAGATTTGCCAATAGGTATTAACATCACT
AGGTTTCAAACCTTACTTGGCTTTACATAGAAGTTATTTGACTCCTGGTGATTCTTCTTCAGGTTGGACAGC
TGGTGTGCAGCTTATTATGTGGGTTATCTTCAACCTAGGACTTTTCTATTAATAATAATGAAAATGGAA
CCATTACAGATGCTGTAGACTGTGCACTTGCACCTCTCTCAGAAAACAAAGTGTACGTTGAAATCCTTCACT
GTAGAAAAAGGAATCTATCAAACCTTCTAAGTTCAGAGTCCAACCAACAGAATCTATTGTTAGATTTCTTAA
TATTACAACTTGTGCCCTTTTGGTGAAGTTTTTAAACGCCACCAGATTTGCATCTGTTTATGCTTGGAA
GGAAGAGAATCAGCAACTGTGTTGCTGATTATTCTGTCTTATATAATTCCGCATCATTTTCCACTTTAAG
TGTTATGGAGTGTCTCTACTAAATTAATGATCTCTGCTTTACTAATGTCTATGCAGATTCATTTGTAAT
TAGAGGTGATGAAGTCAGACAAATCGCTCCAGGGCAAACCTGGAAAGATTGCTGATTATAATTATAAATTAC
CAGATGATTTTACAGGCTGCGTTATAGCTTGGAAATTT**AACA**ATCTTGATTCTAAGGTTGGTGGTAATTTAT
AATTACCTGTATAGATTGTTTAGGAAGTCTAATCTCAAACCTTTTGGAGAGATATTTCAACTGAAATCTA
TCAGGCCGGTAGCACACCTTGAATGGTGTGAAGGTTTTAATTGTTACTTTTCTTTACAATCATATGGTT
TCCAACCCACT**AAT**GGTGTGGTTACCAACCATACAGAGTAGTAGTACTTTCTTTTGAACCTTCTACATGC
ACCAGCAACTGTTTGTGGACCTAAAAAGTCTACTAATTTGGTTAAAAACAAATGTGTCAATTTCAACTTCA
ATGGTTTAAACAGGCACAGGTGTTCTTACTGAGTCTAACAAAAAGTTTCTGCCTTTCCAACAATTTGGCAGA
GACATT**GCT**GACACTACTGATGCTGTCGGTGATCCACAGACACTTGAGATTCTTGACATTACACCATGTT
CTTTTGGTGGTGTGATGTTATAACACCAGGAACAAATACTTCTAACCAGGTTGCTGTTCTTTATCAG**GAT**
GTTAACTGCACAGAAGTCCCTGTTGCTATTCATGCAGATCAACTTACTCCTACTTGGCGTGTTCATTTCTAC
AGGTTCTAATGTTTTTCAAACACGTGCAGGCTGTTAATAGGGGCTGAACATGTCAACAACCTCATATGAGT
GTGACATACCCATTGGTGCAGGTATATGCGCTAGTTATCAGACTCAGACTAATTTCT**CCT**CGGCGGGCAGC
TAGTGTAGCTAGTCAATCCATCATTGCCTACACTATGTCACTTGGTGCAGAAAATTCAGTTGCTTACTCTA
ATAACTCTATTGCCATACCC**ACA**AATTTTACTATTAGTGTACCACAGAAATTTACCAGTGTCTATGAC
CAAGACATCAGTAGATTGTACAATGTACATTTGTGGTGAATCAACTGAATGCAGCAATCTTTTGTGCAAT
ATGGCAGTTTTTGTACACAATTAACCGTGTCTTTAACTGGAATAGCTGTTGAACAAGACAAAAACCCCAA
GAAGTTTTTGCACAAGTCAAACAAATTTACAAAACACCACCAATTAAGATTTTGGTGGTTTTAATTTTTTC
ACAAATATTACCAGATCCATCAAACCAAGCAAGAGGTCATTTATTGAAGATCTACTTTTCAACAAGTGA
CACTTGCAGATGCTGGCTTCATCAAACAATATGGTGAATGCCTTGGTGAATGCTGCTAGAGACCTCATT
TGTGCACAAAAGTTTTAACGGCCTTACTGTTTTGCCACCTTTGCTCACAGATGAAATGATTGCTCAATACAC
TTCTGCACTGTTAGCGGGTACAATCACTTCTGGTTGGACCTTTGGTGCAGGTGCTGCATTACAAATACCAT
TTGCTATGCAAATGGCTTATAGGTTAATGGTATTGGAGTTACACAGAATGTTCTCTATGAGAACCAAAAA
TTGATTGCCAACCAATTTAATAGTGTCTATTGGCAAAATTCAGACTCACTTTCTTCCACAGCAAGTCACT
TGGAAAACCTCAAGATGTGGTCAACCAAAATGCACAAGCTTTAAACACGCTTGTAAACAACCTTAGCTCCA
ATTTTGGTGAATTTCAAGTGTTTTAAATGATATCCTT**TCA**CGTCTTGACAAAGTTGAGGCTGAAGTGCA
AATTGATAGGTTGATCACAGGCAGACTTCAAAGTTTGCAGACATATGTGACTCAACAATTAATTAGAGCTG
CAGAAATCAGAGCTTCTGCTAATCTTGCTGCTACTAAAATGTGAGAGTGTGACTTGGACAATCAAAAA
GTTGATTTTTTGTGAAAGGGCTATCATCTTATGTCTTCCCTCAGTCAGCACCTCATGGTGTAGTCTTCTT
GCATGTGACTTATGTCCCTGCACAAGAAAAGAACTTCAACTGCTCCTGCCATTTGTCATGATGGAAAAG
CACACTTTCTCGTGAAGGTGTCTTTGTTTCAAATGGCACACACTGGTTTGTAAACACAAAGGAATTTTAT
GAACCACAAATCATTACTACA**GAC**AACACATTTGTGTCTGGTAACTGTGATGTTGTAATAGGAATTGTCA
ACAACACAGTTTATGATCCTTTGCAACCTGAATTAGACTCATTCAAGGAGGAGTTAGATAAATATTTAAG
AATCATACATCACCAGATGTTGATTTAGGTGACATCTCTGGCATTAAATGCTTCAAGTTGTAACATTCAAAA
AGAAATTGACCGCCTCAATGAGGTTGCCAAGAATTTAAATGAATCTCTCATCGATCTCCAAGAACTTGAA
AGTATGAGCAGTATATAAATGGCCATGGTACATTTGGCTAGGTTTTATAGCTGGCTTGATTGCCATAGTA

ATGGTGACAATTATGCTTTGCTGTATGACCAGTTGCTGTAGTTGTCTCAAGGGCTGTTGTTCTTGTGGATC
CTGCTGCAAATTTGATGAAGACGACTCTGAGCCAGTGCTCAAAGGAGTCAAATTACATTACACATAA

3. Conclusions

After the analysis of the sequence, GeneMe can state that:

- The observed mutations within the S gene are outside the sequence of primers utilized in the FRANKD test. Sequences of the primers are a trade secret of GeneMe; therefore, they are not disclosed in this report. Upon specific non-disclosure agreements, GeneMe can reveal those sequences to the regulatory bodies – notified bodies.

Further notice:

- It is important to remember that point mutations even in the hybridization region of primer sequences may not significantly influence the efficacy of the primer's hybridization.
- FRANKD, as an RT-LAMP test uses six primers for identification rather than the two normally used by RT-PCR tests.
- FRANKD utilized its own patented Phusion polymerase with the SSB domain that guarantees stronger binding of the primers to the target sequence after the hybridization, limiting the possibility of temperature influence on the dissociation of the primers.
- GeneMe constantly follows scientific information on point mutations in the S gene. Upon any new information, the Head of Research and Development performs *in silico* analysis of the hybridization efficacy of FRANKD primers to the S gene of SARS-CoV-2. In case of a point mutation or a deletion that may influence the test's performance, a new primer set will be designed by changing the point mutated nucleotide for a degenerated nucleotide.
- GeneMe, as a company, manages five diagnostic laboratories in Poland (laboratories certified by the Polish Health Ministry). As a service, GeneMe provides patient sample analysis for the presence of SARS-CoV-2 with the use of RT-PCR and RT-LAMP kits targeting E, Orf1ab, E and S genes. Therefore, GeneMe has the possibility to monitor the

performance of the diagnostic kits available on the market and produced by GeneMe. It was already established that point mutations also occur in other target sequences as presented below [5, 6, 7]:

| gene | nucleotide | amino acid | |
|---------|----------------------|------------------------|---------|
| ORF1ab | C3267T | T1001I | |
| | C5388A | A1708D | |
| | T6954C | I2230T | |
| | 11288-11296 deletion | SGF 3675-3677 deletion | |
| spike | 21765-21770 deletion | HV 69-70 deletion | |
| | 21991-21993 deletion | Y144 deletion | |
| | A23063T | N501Y | |
| | C23271A | A570D | |
| | C23604A | P681H | |
| | C23709T | T716I | |
| | T24506G | S982A | |
| | G24914C | D1118H | |
| | Orf8 | C27972T | Q27stop |
| | | G28048T | R52I |
| A28111G | | Y73C | |
| N | 28280 GAT->CTA | D3L | |
| | C28977T | S235F | |

- At GeneMe, we are aware of the possibility of fast mutations that can occur in the group of coronaviruses. Therefore GeneMe has established procedures to monitor those changes and their potential influence on the performance of diagnostic kits produced by us.
- **In our opinion, RT-LAMP technologies are less prone to lose specificity and sensitivity in comparison to RT-PCR, because of (i) amplification principle and reaction efficacy (constant amplification, no doubling DNA in each cycle, (ii) more than two primers in the reaction (more targets), (iii) polymerase features (strand display activity of the polymerase) [8].**

Reference to the literature:

- 
- [1] https://www.cogconsortium.uk/news_item/update-on-new-sars-cov-2-variant-and-how-cog-uk-tracks-emerging-mutations/
- [2] Thomson EC, Rosen LE, Shepherd JG, et al. The circulating SARS-CoV-2 spike variant N439K maintains fitness while evading antibody-mediated immunity. *bioRxiv* 2020: 2020.11.04.355842.
- [3] Weissman D, Alameh M-G, de Silva T, et al. D614G Spike Mutation Increases SARS CoV-2 Susceptibility to Neutralization. *medRxiv* 2020: 2020.07.22.20159905.
- [4] [Rapid increase of a SARS-CoV-2 variant with multiple spike protein mutations observed in the United Kingdom](#)
- [5] [Preliminary genomic characterization of an emergent SARS-CoV-2 lineage in the UK defined by a novel set of spike mutations](#)
- [6] Khailany RA, Safdar M, Ozaslan M. Genomic characterization of a novel SARS-CoV-2. *Gene Rep.* 2020;19:100682. doi:10.1016/j.genrep.2020.100682
- [7] [A Recurrent Mutation at Position 26340 of SARS-CoV-2 Is Associated with Failure of the E Gene Quantitative Reverse Transcription-PCR Utilized in a Commercial Dual-Target Diagnostic Assay](#)
- Maria Artesi, Sébastien Bontems, Paul Göbbels, Marc Franckh, Piet Maes, Raphaël Boreux, Cécile Meex, Pierrette Melin, Marie-Pierre Hayette, Vincent Bours, Keith Durkin
- Journal of Clinical Microbiology* Sep 2020, 58 (10) e01598-20; DOI: 10.1128/JCM.01598-20
- [8] [Olszewski M, Balsewicz J, Nowak M, Maciejewska N, Cyranka-Czaja A, Zalewska-Piątek B, et al. \(2015\) Characterization of a Single-Stranded DNA-Binding-Like Protein from *Nanoarchaeum equitans*—A Nucleic Acid Binding Protein with Broad Substrate Specificity. *PLoS ONE* 10\(5\): e0126563. <https://doi.org/10.1371/journal.pone.0126563>](#)

Prepared

RD Head [Marta Skwarecka](#)

M. Skwarecka

Authorized

RD Director [Kasjan Szemiako](#)

K. Szemiako